

Chifa

Protetyka i Ortodoncja

Prosthetics and Orthodontics



New models

Torquing pliers and plastic ligature applying tool

Excellence in OEM manufacturing


www.chifa-oem.pl




Kleszcze do torkowania łuków Torquing pliers

KO-914-140-PZMK

140 mm
max 0,55 x 0,64 mm
szerokość torku 6,5 mm
arc width 6,5mm (0.0256")




1/1



Torquing scale on working end inform of the arc scale

KO-915-140-PZMK

140 mm
max 0,55 x 0,64 mm
szerokość torku 6,5 mm
arc width 6,5mm (0.0256")



1/1



Wider working part for easier fit between brackets
and to create the arc between brackets in a single move

Narzędzie do zaginania łuków
Instrument for bending arch wire ends

SP-009-155-PMK

155 mm
max -0,8 mm



Instrument for bending (90°) arch wire ends

Sierp ortodontyczny z pozycjonerem
Instrument for applying plastic ligature

SP-014-160-PMK

160 mm



With measuring function

