



Basing on our international experience we designed a unique line of aligners. The pliers meet the expectations of orthodontists providing optimal support of the work. Made from high grade stainless steel along with precise workmanship, guarantee perfect functionality.

CHIFA[®]

NEW LINE OF ALIGNER PLIER




KO-125-120-PMK

 120 mm

The Horizontal.

Plier enables to horizontally correct clear aligner tooth positioners with no punching




KO-126-120-PMK


 120 mm

The Vertical.

Plier enables to vertically correct clear aligner tooth positioners with no punching




KO-128-120-PMK

 120 mm

The Hole Punch

The plier produces half-moon shape to allow attachments to be placed.



KO-129-120-PMK

 120 mm

The Tear Drop.

Creates a notch in clear aligner in an area, where hooks are required.

Aesculap-Chifa Sp. z o.o. Tysiąclecia 14. 64-300-Nowy Tomyśl, Poland

tel. +48 61 44 20 261, fax. 61 44 20260

info@chifa-oem.pl www.chifa-oem.pl